

Rad-57™

OPERATOR'S MANUAL

Masimo Rainbow® SET®

Signal Extraction
Pulse CO-Oximeter™

- Covers:
 - Rad-57c:
with SpO₂ and SpCO
 - Rad-57m:
with SpO₂ and SpMet
 - Rad-57cm:
with SpO₂, SpCO
and SpMet



Masimo SET®
rainbow®

The Rad-57 Signal Extraction Pulse CO-Oximeter Operating Instructions intend to provide the necessary information for proper operation of all Rad-57 Pulse CO-Oximeter models.

General knowledge of Pulse CO-Oximetry and an understanding of the features and functions of the Rad-57 Signal Extraction Pulse CO-Oximeter models are prerequisites for proper use.

Do not operate any of the Rad-57 Signal Extraction Pulse CO-Oximeter models without completely reading and understanding these instructions.

NOTICE

Purchase or possession of this device does not carry any express or implied license to use with replacement parts which would, alone or in combination with this device, fall within the scope of one of the relating patents.

CAUTION:

FEDERAL LAW (U.S.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For further information contact:

Masimo Corporation

40 Parker

Irvine, CA 92618

USA


Tel.: 949-297-7000

Fax.: 949-297-7001

www.masimo.com



EU Authorized Representative for Masimo Corporation:



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany
Tel.: +49-511-62 62 86 30
Fax.: +49-511-62 62 86 33

Covered by one or more of the following U.S. Patents: RE38,492, RE38,476, 6,699,194, 6,684,090, 6,658,276, 6,654,624, 6,650,917, 6,643,530, 6,606,511, 6,501,975, 6,463,311, 6,430,525, 6,360,114, 6,263,222, 6,236,872, 6,229,856, 6,206,830, 6,157,850, 6,067,462, 6,011,986, 6,002,952, 5,919,134, 5,823,950, 5,769,785, 5,758,644, 5,685,299, 5,632,272, 5,490,505, 5,482,036, international equivalents, or one or more of the patents referenced at www.masimo.com/patents. Other patents pending.

Manufactured in USA

© 2006 Masimo Corporation. Masimo, SET, LNCS, LNOP, Rainbow, SpCO, Signal IQ, Discrete Saturation Transform, DST and FasiSat are registered trademarks of Masimo Corporation. Rad-57, Rad-57c, Rad-57m, Rad-57cm, SIQ, LNOPv, APOD, SpMet and Pulse CO-Oximeter are trademarks of Masimo Corporation.

SAFETY INFORMATION, WARNINGS, CAUTIONS AND NOTES

The Rad-57™ Handheld Pulse CO-Oximeter™ is designed to minimize the possibility of hazards from errors in the software program by following sound engineering design processes, Risk Analysis and Software Validation.

- Explosion hazard. Do not use the Pulse CO-Oximeter in the presence of flammable anesthetics or other flammable substance in combination with air, oxygen-enriched environments, or nitrous oxide.
- High intensity extreme lights (such as pulsating strobe lights) directed on the sensor, may not allow the Pulse CO-Oximeter to obtain vital sign readings.
- The Pulse CO-Oximeter is NOT intended for use as an apnea monitor.
- A Pulse CO-Oximeter should be considered an early warning device. As a trend towards patient hypoxemia is indicated, blood samples should be analyzed by laboratory instruments to completely understand the patient's condition.
- The Pulse CO-Oximeter is to be operated by qualified personnel only. This manual, accessory directions for use, all precautionary information, and specifications should be read before use.
- Electric shock hazard. Do not open the Pulse CO-Oximeter cover except to replace the batteries. Only a qualified operator may perform maintenance procedures specifically described in this manual. Refer servicing to Masimo for repair of this equipment.
- As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.
- Do not place the Pulse CO-Oximeter or accessories in any position that might cause it to fall on the patient. Do not lift the Pulse CO-Oximeter by the patient cable.
- Interfering Substances: SpO₂ is a functional calculation of arterial oxygen saturation. Carboxyhemoglobin and methemoglobin may erroneously increase SpO₂ readings. The level of increase is approximately equal to the amount of carboxyhemoglobin and/or methemoglobin present. Dyes, or any substance containing dyes, that change usual blood pigmentation may cause erroneous readings.
- Severe anemia may cause erroneous SpO₂ readings.
- Do not use the Pulse CO-Oximeter or oximetry sensors during magnetic resonance imaging (MRI) scanning. Induced current could potentially cause burns. The Pulse CO-Oximeter may affect the MRI image, and the MRI unit may affect the accuracy of the oximetry measurements.
- If using Pulse CO-Oximetry during full body irradiation, keep the sensor out of the irradiation field. If sensor is exposed to the irradiation, the reading might be inaccurate or the unit might read zero for the duration of the active irradiation period.
- Always remove the sensor from the patient and completely disconnect the patient from the Pulse CO-Oximeter before bathing the patient.
- Do not place the Pulse CO-Oximeter where the controls can be changed by the patient.
- Do not place the Pulse CO-Oximeter face against a surface. This will cause the alarm to be muffled.
- Do not place the Pulse CO-Oximeter on electrical equipment that may affect the Pulse CO-Oximeter, preventing it from working properly.

SAFETY INFORMATION, WARNINGS, CAUTIONS AND NOTES

- Do not expose the Pulse CO-Oximeter to excessive moisture such as direct exposure to rain. Excessive moisture can cause the Pulse CO-Oximeter to perform inaccurately or fail.
- Do not place containers containing liquids on or near the Pulse CO-Oximeter. Liquids spilled on the Pulse CO-Oximeter may cause it to perform inaccurately or fail.
- Failure of Operation - If the Pulse CO-Oximeter fails any part of the setup procedures remove the Pulse CO-Oximeter from operation until qualified service personnel have corrected the situation.
- Patient Safety - If a sensor is damaged in any way, discontinue use immediately.
- Disposal of product - Comply with local laws in the disposal of the unit and/or its accessories.
- The Pulse CO-Oximeter can be used during defibrillation, but the readings may be inaccurate for up to 20 seconds.
- This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2: 2002, Medical Device Directive 93/42/EEC. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving device.
 - Increase the separation between the equipment.
 - Consult the manufacturer for help.
- A functional tester cannot be utilized to assess the accuracy of the Pulse CO-Oximeter or any sensors.

table of contents

SECTION 1 - OVERVIEW	1-1
About this Manual.....	1-1
Warnings, Cautions and Notes.....	1-2
Product Description	1-3
Model Summary	1-3
Features and Benefits	1-3
Indications for Use	1-4
Pulse CO-Oximetry	1-5
SpO ₂ General Description	1-5
SpCO General Description.....	1-5
SpMet General Description	1-5
Principle of Operation	1-6
Measured vs. Calculated Values	1-7
Masimo SET Signal Extraction Technology for SpO ₂ measurements.....	1-7
SpCO and SpMet Measurements During Patient Motion	1-7
Masimo SET Parallel Engines.....	1-8
Masimo SET DST®	1-8
SECTION 2 - SYSTEM DESCRIPTION	2-1
Introduction.....	2-1
Rad-57c Front Panel Controls	2-2
Rad-57m Front Panel Controls	2-4
Rad-57cm Front Panel Controls	2-6
Symbols.....	2-8
Rad-57 Rear Panel (all models).....	2-8
SECTION 3 - SETUP	3-1
Introduction.....	3-1
Preparation for monitoring.....	3-1
Power requirements	3-1
Monitor setup.....	3-2
Initial Setup.....	3-2
SECTION 4 - OPERATION	4-1
Introduction.....	4-1
Basic operation.....	4-1
General Setup and Use.....	4-1
Default settings.....	4-3
Successful SpO ₂ monitoring	4-4
Numeric Display - SpO ₂	4-4
Masimo Sensors.....	4-4
Numeric Display - Pulse Rate	4-5
Numeric Display - SpCO (for applicable models).....	4-5
Numeric Display - SpMet (for applicable models).....	4-5
Low Signal IQ (Low SIO).....	4-5
Low Perfusion.....	4-6
Actions To Be Taken.....	4-6
Sensitivity	4-6
Low Battery Audible Alarm.....	4-8

table of contents

Normal Patient Monitoring.....	4-9
Setup Menu.....	4-10
Menu Navigation	4-10
Setup Menu Level 2 - Alarm Limits	4-10
Setup Menu Level 3 - Averaging and Sensitivity	4-11
Setup menu Level 4 - Trend Settings	4-11
Setup Menu Level 5 - LED Brightness and Factory Defaults.....	4-12
Menu Selection.....	4-12
Power Off.....	4-12
Special Menu.....	4-13
Special Menu - Line Frequency Configuration.....	4-13
Trend Setup and Use.....	4-14
Introduction.....	4-14
TrendCom Utility Installation.....	4-14
TrendCom Utility operation.....	4-14
Erasing Trend Memory	4-15
Trend Data Format	4-15
SECTION 5 - ALARMS / MESSAGES	5-1
Alarm Indication.....	5-1
Alarm Limits.....	5-1
Alarm Silence	5-3
Alarm Silenced Indicator	5-3
SECTION 6 - TROUBLESHOOTING	6-1
Troubleshooting	6-1
SECTION 7 - SPECIFICATIONS	7-1
Rad-57 Family Specifications	7-1
Performance.....	7-1
Accuracy.....	7-1
Electrical.....	7-2
Environmental.....	7-2
Physical characteristics	7-2
SECTION 8 - SENSOR AND PATIENT CABLES	8-1
Introduction.....	8-1
Selecting a Sensor	8-1
Sensor Application Site	8-1
Masimo Rainbow® Sensors	8-2
Rainbow Reusable Sensors	8-2
Masimo SpO ₂ Sensors.....	8-2
Red Reusable Sensors	8-2
LNOP® Reusable Sensors	8-2
LNOP® Adhesive Sensors	8-3
LNOP™ Specialty Sensors	8-3
LNCS™ Reusable Sensors	8-3
LNCS™ Adhesive Sensors	8-3
LNOPV™ Adhesive Sensors	8-4
Cleaning And Reuse of Masimo Reusable Sensors and Cables.....	8-4
Reattachment of Single Use Adhesive Sensors	8-4

SECTION 9 - SERVICE / MAINTENANCE..... 9-1

Introduction..... 9-1

Cleaning..... 9-1

Battery Replacement..... 9-2

Performance Verification..... 9-3

Power-On Self-Test..... 9-3

Key Press Button Test..... 9-3

Alarm Limit Test..... 9-3

LED Brightness..... 9-4

Service and Repair..... 9-4

 Repair Policy..... 9-4

 Return Procedure..... 9-4

Warranty..... 9-5

Exclusions..... 9-5

End-user License Agreement..... 9-6

SECTION 10 - ACCESSORIES..... 10-1

Accessories..... 10-1

About this Manual

This manual explains how to set up and use the Rad-57c, Rad-57m and Rad-57cm Pulse CO-Oximeter. These three products will be generally referred to as the Rad-57. Features and specifications apply to all models unless otherwise noted. Important safety information relating to general use of the Rad-57 appears before this introduction. Other important safety information is located throughout the manual where appropriate.

Read the entire safety information section before you operate the monitor.

In addition to the safety section, this manual includes the following sections:

- SECTION 1** OVERVIEW gives a general description of pulse oximetry.
- SECTION 2** SYSTEM DESCRIPTION describes the Rad-57 system, functions and features.
- SECTION 3** SETUP describes how to setup the Rad-57 for use.
- SECTION 4** OPERATION describes the operation of the Rad-57.
- SECTION 5** ALARMS AND MESSAGES describes the alarm system messages.
- SECTION 6** TROUBLESHOOTING gives troubleshooting information.
- SECTION 7** SPECIFICATIONS gives the detailed specifications of the Rad-57.
- SECTION 8** SENSORS AND PATIENT CABLES outlines how to use and care for the Rainbow, LNOP®, LNOPY™ and LNCS® sensors and Masimo SET patient cables.
- SECTION 9** SERVICE AND MAINTENANCE describes how to maintain, service and obtain repair for the Rad-57.
- SECTION 10** ACCESSORIES

Warnings, Cautions and Notes

Please read and follow any warnings, cautions and notes presented throughout this manual. An explanation of these labels are as follows:

A **WARNING** is provided when actions may result in a serious outcome (i.e., injury, serious adverse affect, death) to the patient or user. Look for text in a gray shaded box. Sample of Warning:

WARNING: THIS IS A SAMPLE OF A WARNING STATEMENT.

A **CAUTION** is given when any special care is to be exercised by the patient or user to avoid injury to the patient, damage to this device or damage to other property. Sample of Caution:

CAUTION: THIS IS A SAMPLE OF A CAUTION STATEMENT.

A **NOTE** is provided when extra general information is applicable. Sample of Note:

NOTE: *This is a sample of a Note.*

Product Description

The Rad-57 Handheld Pulse CO-Oximeter with Masimo Rainbow® SET® Technology is a noninvasive, arterial oxygen saturation and pulse rate monitor. The Rad-57 features a multicolored LED display that continuously displays numeric values for SpO₂ and pulse rate, a Low Signal IQ Indicator (Low SIQ) indicator, alarm status, alarm silence and battery life.

Certain models include LED indicator bars for Perfusion Index (PI), Carboxyhemoglobin saturation (%SpCO), Methemoglobin percentage (%SpMet), and SpCO/SpMet sensor connected.

MODEL SUMMARY

All models include SpO₂ and pulse rate.

MODEL	FEATURES
Rad-57c	Includes SpCO and PI
Rad-57m	Includes SpMet and PI
Rad-57cm	Includes SpCO, SpMet and Low PI Indicator

The following list outlines the key features and benefits of the Rad-57 Handheld Pulse CO-Oximeter.

FEATURES AND BENEFITS

- Clinically proven Masimo SET® technology performance
- Proven for accurate monitoring in motion and low perfusion environments
- SpO₂, pulse rate, alarm, Perfusion Index*, % SpCO* and SpMet* displays
- Low Signal IQ (SIQ®) indicator
- Lightweight, convenient handheld design
- Up to 8 hours of continuous use on 4 "AA" alkaline batteries
- Visual battery life indicator
- Audible Alarm for sensor-off and low battery
- Alarms for high/low SpO₂, high/low pulse rate, SpCO* and SpMet*
- FastSat® (for SpO₂ measurement)
- Three sensitivity levels - Max, Normal and APOD™ (for SpO₂ measurement)
- 72 hours of trending memory
- Adjustable alarm volume
- Adjustable averaging 2 to 16 seconds

*See Model Summary for applicable unit.

INDICATIONS FOR USE

The Rad-57 and accessories are indicated for the continuous, noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂), pulse rate, carboxyhemoglobin saturation (measured by an SpCO/SpMet sensor)* and methemoglobin percentage (measured by a SpCO/SpMet sensor)*. The Rad-57 and accessories are indicated for use with adult, pediatric and neonatal patients during both motion and no motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile and home environments.

*See Model Summary for applicable unit.


Pulse CO-Oximetry

SpO₂ GENERAL DESCRIPTION

Pulse oximetry is a continuous and non-invasive method of measuring the level of arterial oxygen saturation in blood. The measurement is taken by placing a sensor on a patient, usually on the fingertip for adults, and the hand or foot for neonates. The sensor connects to the pulse oximetry instrument with a patient cable. The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data in two ways:

- 1) As a percent value for arterial oxygen saturation (SpO₂)
- 2) As a pulse rate (PR)

The following figure shows the general monitoring setup.



SpCO GENERAL DESCRIPTION

Pulse CO-Oximetry is a continuous and non-invasive method of measuring the levels of carbon monoxide concentration (SpCO) in arterial blood. It relies on the same principles of pulse oximetry to make its SpCO measurement. The measurement is taken by placing a sensor on a patient, usually on the fingertip for adults. The sensor connects directly to the pulse CO-Oximetry instrument or with a patient cable. The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data as percentage value for the SpCO. The Rad-57c is a combined SpO₂ and SpCO monitor with the same setup as that of a pulse oximeter, shown above, and can display a percentage value for SpCO as well as SpO₂. The Rad-57cm displays the same parameters as the Rad-57c and also displays SpMet.


SpMet GENERAL DESCRIPTION

Pulse CO-Oximetry is a continuous and non-invasive method of measuring the levels of methemoglobin concentration (SpMet) in arterial blood. It relies on the same principles of pulse oximetry to make its SpMet measurement. The measurement is taken by placing a sensor on a patient, usually on the fingertip for adults. The sensor connects directly to the Pulse CO-Oximetry instrument or with a patient cable. The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data as percentage value for the SpMet. The Rad-57m is a combined SpO₂ and SpMet monitor with the same setup as that of a pulse oximeter, shown above, and can display a percentage value for SpMet as well as SpO₂. The Rad-57cm displays the same parameters as the Rad-57m and also displays SpCO.

PRINCIPLE OF OPERATION


Pulse oximetry is governed by the following principles:

1. Oxyhemoglobin (oxygenated blood), deoxyhemoglobin (non-oxygenated blood), carboxyhemoglobin (blood with carbon monoxide content) and methemoglobin (blood with oxidized hemoglobin content) species differ in their absorption of visible and infrared light (using spectrophotometry, see figure below).



2. The amount of arterial blood in tissue changes with your pulse (photoplethysmography). Therefore, the amount of light absorbed by the varying quantities of arterial blood changes as well.

The Rad-57 Pulse CO-Oximeter uses a multi-wavelength sensor to distinguish between oxygenated blood, deoxygenated blood, blood with carbon monoxide content and blood with oxidized hemoglobin content. Signal data is obtained by passing various visible and infrared lights (LEDs, 500nm to 1000nm) through a capillary bed (for example, a fingertip, a hand, a foot) and measuring changes in light absorption during the blood pulsatile cycle. This information may be useful to clinicians. The maximum radiant power of the strongest light is rated at 22mW. See figure below. The photodetector receives the light, converts it into an electronic signal and sends it to the Rad-57 for calculation.



Once the Rad-57 receives the signal from the sensor, it utilizes Masimo Rainbow SET Technology signal extraction technology to calculate the patient's functional oxygen saturation, fractional concentration of carboxyhemoglobin, methemoglobin and pulse rate. The SpCO and SpMet measurements rely on a multiwavelength calibration equation to estimate the percentage of carbon monoxide and oxidized hemoglobin in arterial blood. The maximum of the skin surface temperature is measured at an ambient temperature of less than 106° F (41° C). This is verified by Masimo sensor skin temperature test procedures.

FUNCTIONAL VS. FRACTIONAL SATURATION

The Rad-57 is calibrated to measure and display *functional saturation* (SpO₂): the amount of oxyhemoglobin expressed as a percentage of the hemoglobin that can transport oxygen. The Rad-57 does not measure *fractional saturation*: oxygenated hemoglobin expressed as

a percentage of the four main hemoglobin species, i.e., oxyhemoglobin, deoxyhemoglobin, carboxyhemoglobin, and methemoglobin. To convert fractional saturation to functional saturation, the fractional saturation measurements must be converted according to:

$$\text{Functional saturation} = \frac{\text{Fractional saturation}}{100 - (\% \text{ carboxyhemoglobin} + \% \text{ methemoglobin})} \times 100$$

MEASURED VS. CALCULATED VALUES

SpO₂, SpCO and SpMet measurements can be obtained from the Rad-57. Each functional parameter corresponds to the applicable model. They are commonly compared to invasive measurements obtained from blood gas samples. When comparing invasive and noninvasive measurements and interpreting values, caution should be used, as the calculated values obtained from the blood gas sample may differ from the SpO₂, SpCO and SpMet measurements of the Pulse CO-Oximeter. In the case of SpO₂, different results are usually obtained from the arterial blood gas sample if the calculated measurement is not appropriately corrected for the effects of variables that shift the relationship between the partial pressure of oxygen (PO₂) and saturation, such as: pH, temperature, the partial pressure of carbon dioxide (PCO₂), 2,3-DPG, and fetal hemoglobin. In the case of SpCO and SpMet, in addition to the effects of temperature and pH, different results are also expected if the oxygen saturation and/or concentration of methemoglobin in the blood gas sample are abnormal (less than 90% for arterial oxygen saturation, and greater than 1% for methemoglobin concentration). As blood gas samples are usually taken over a period of 20 seconds (the time it takes to draw the blood) a meaningful comparison can only be achieved if the core oxygen saturation, methemoglobin concentration and carboxyhemoglobin concentration of the patient are stable and not changing over the period of time that the blood gas sample is taken.

MASIMO SET SIGNAL EXTRACTION TECHNOLOGY FOR SpO₂ MEASUREMENTS

Masimo Signal Extraction Technology's signal processing differs from conventional pulse oximeters. Conventional pulse oximeters assume that arterial blood is the only blood moving (pulsating) in the measurement site. During patient motion, however, the non-arterial blood also moves, causing conventional pulse oximeters to read low values, because they cannot distinguish between the arterial and venous blood movement (sometimes referred to as noise). Masimo SET pulse oximetry utilizes parallel engines and adaptive digital filtering. Adaptive filters are powerful because they are able to adapt to the varying physiologic signals and/or noise and separate them by looking at the whole signal and breaking it down to its fundamental components. The Masimo SET signal processing algorithm, Discrete Saturation Transform® (DST®), reliably identifies the noise, isolates it and, using adaptive filters, cancels it. It then reports the true arterial oxygen saturation for display on the monitor.

SpCO AND SpMet MEASUREMENTS DURING PATIENT MOTION

Rad-57 displays measurements of SpCO and SpMet (each parameter corresponding to the applicable model) during patient motion. However, because of the changes in the physiological parameters such as blood volume, arterial-venous coupling, etc. that occur during patient motion, the accuracy of such measurements are not reliable.

1

overview

MASIMO SET PARALLEL ENGINES

The flowchart illustrates the parallel processing pipeline for MASIMO SET. It begins with an input of $R/I/R$ which is processed by a 'Digitized, Filtered & Normalized' block. This signal is then fed into five parallel processing engines: 'R/I/R (Conventional Pulse Oximetry)', 'Adaptive Filter with DST[®]', 'SST[™]', 'Proprietary Algorithm 4', and 'Proprietary Algorithm 5'. Each engine outputs a 'MEAS' (Measurement) and a 'CONF' (Confidence) signal. These signals are combined in a 'Confidence Based Arbitrator' block. The output of the arbitrator goes to a 'Post Processor', which finally outputs the 'Saturation' level.

MASIMO SET DST[®]

The diagram shows the MASIMO SET DST process. It starts with 'IR' and 'RD' signals, which are processed by a 'Reference Signal Generator'. This generator also receives 'Trial Saturation SpO₂ 95%' as input. The output is a 'Reference Signal' that is fed into an 'Adaptive Filter'. The filter's output is 'Output Power', which is plotted on a graph. The graph shows a bell-shaped curve of 'Power' versus 'SpO₂ (r_v)' on the x-axis (ranging from 30 to 100) and 'SpO₂ (r_a)' on the y-axis. The peak of the curve is at approximately 95% saturation.

system description

2

Introduction

The Rad-57 is a full featured Pulse CO-Oximeter designed for ease of operation. All pulse oximetry measurement information, as well as device status data, is displayed on the front panel of the device. All user input is handled by control buttons on the front panel and the sensor cable connection is located at the top edge of the device.

The Rad-57 is powered by 4 "AA" alkaline batteries, which provides a minimum of 8 hours of battery life.

- Rad-57 offers full Masimo SET with Rainbow Technology in a small, hand held device.
- Rad-57 supports the full line of Masimo sensors (see Section 8, *sensors and patient cables*).
- Provides 72 hours of trending memory

A Rainbow DC-I/DC-IP Sensor or a Red Patient Cable with a LNCS, LNCP or LNCPv sensor attaches to the patient cable connector on the top of the Rad-57 unit. The Rad-57 can be used either as a transport monitor or as a handheld Pulse CO-Oximeter for spot checks.


1-8

Rad-57 Signal Extraction CO-Pulse Oximeter Operator's Manual

Rad-57 Signal Extraction CO-Pulse Oximeter Operator's Manual

2-1

Rad-57c Front Panel Controls




Rad-57c

CONTROL / INDICATOR	DESCRIPTION
① Patient Cable Connector	Connects to a Rainbow sensor, Red sensor or a Red Patient Cable with an LNOP, LNOPv or LNCS sensor.
② SpCO Indicator	Slow flashing indicator: The confidence in the SpCO value obtained is low. Fast flashing indicator: Flashes when an SpCO alarm condition exists.
③ Saturation (%SpO ₂) and Carboxyhemoglobin (%SpCO) Displays	The functional arterial hemoglobin oxygen saturation is displayed in units of SpO ₂ . When searching for a saturation and pulse, it will flash dashed lines. See SpCO Button description, below.
④ %SpCO Bar	Illuminates when SpCO capable sensor is attached. Bar will flash for SpCO alarm conditions. Continuously indicates the concentration of carboxyhemoglobin in 5% increments.
⑤ Pulse Rate Display	The pulse rate in beats per minute (bpm). When searching for a saturation and pulse, it will flash dashed lines.

CONTROL / INDICATOR	DESCRIPTION
⑥ SpCO Button	Pressing this button will display the numeric SpCO value for 10 seconds in place of the SpO ₂ numeric value. Pressing the Mode/Enter or Next button during this 10-second period will return to the SpO ₂ numeric value.
⑦ Mode / Enter Button	Used to enter the setup menus and to select/activate certain entries within the menu/setup system.
⑧ Next Button	Used within the menu/setup system to move through setup options. Not active during normal patient monitoring.
⑨ Power On / Off	Press to turn the unit on. Press-and-hold for 2 seconds to turn the unit off.
⑩ Battery Level Indicator	Four LED's indicate the status of the battery. When the final indicator begins flashing, replace the batteries.
⑪ Visual Alarm Indicator	Illuminates when any alarm condition exists. This indicator may not be turned-off or otherwise over-ridden.
⑫ %SpO ₂ / %SpCO Indicator	Indicator above label will illuminate to provide an additional visual indication of the value currently being displayed.
⑬ PI	Perfusion Index, or PI, is a relative assessment of the perfusion at the monitoring site. PI is displayed on a 10 segment LED bar, ranging from <1% (very weak perfusion) to >5% (strong perfusion). The highest LED will remain lit continuously to allow a PI level to be viewed. The Perfusion Index is the ratio of the AC (pulsatile) to DC (non-pulsatile) components of the IR (Infrared) signal where the AC and DC components correspond to the pulsatile and non-pulsatile amounts of blood, respectively.
⑭ Low SIQ	Flashes to indicate low SpO ₂ Signal IQ. Refer to Section 4, Low Signal IQ.
⑮ Alarm Silenced Indicator	Flashes to indicate the alarm is temporarily silenced.
⑯ Alarm Silence Button	Push once to temporarily silence the alarm for 120 seconds. Push a second time to return the unit to standard alarm monitoring.
⑰ Up button Down button	During saturation monitoring, use these buttons to adjust the volume of the pulse beep tone. Within the menu/setup system, these buttons are used to select values within each menu option.
⑱ Speaker	Provides audible indication of alarm conditions, pulse tone and feedback for key-presses. Ensure the speaker is not covered or the unit is placed face-down on bedding or other sound absorbing surface.

Rad-57m Front Panel Controls




Rad-57m

CONTROL / INDICATOR	DESCRIPTION
① Patient Cable Connector	Connects to a Rainbow sensor, Red sensor or a Red Patient Cable with an LNOP, LNOPv or LNCS sensor.
② SpMet indicator	Slow flashing indicator: The confidence in the SpMet value obtained is low. Fast flashing indicator: Flashes when an SpMet alarm condition exists.
③ Saturation (%SpO ₂) Display	The functional arterial hemoglobin oxygen saturation is displayed in units of SpO ₂ . When searching for a saturation and pulse, it will flash dashed lines.
④ %SpMet Bar	Illuminates when SpMet capable sensor is attached. Bar will flash for SpMet alarm conditions. Continuously indicates the concentration of Methemoglobin in increments: .5, 1-5, 7-5, 10, 15 >20%.
⑤ %SpMet / Pulse Rate Display	Indicator above label will illuminate to provide an additional visual indication of the SpMet value currently being displayed. The pulse rate in beats per minute (bpm). When searching for a saturation and pulse, it will flash dashed lines.

CONTROL / INDICATOR	DESCRIPTION
⑥ SpMet Button	Pressing this button will display the numeric SpMet value in place of the SpO ₂ numeric value until the SpMet button is pressed again.
⑦ Mode / Enter Button	Used to enter the setup menu and to select/activate certain entries within the menu/setup system.
⑧ Next Button	Used within the menu/setup system to move through setup options. Not active during normal patient monitoring.
⑨ Power On / Off	Press to turn the unit on. Press-and-hold for 2 seconds to turn the unit off.
⑩ Battery Level Indicator	Four LED's indicate the status of the battery. When the final indicator begins flashing, replace the batteries.
⑪ Visual Alarm Indicator	Illuminates when any alarm condition exists. This indicator may not be turned-off or otherwise over-ridden.
⑫ %SpO ₂ Indicator	Indicator above label will illuminate to provide an additional visual indication of the value currently being displayed.
⑬ PI	Perfusion Index, or PI, is a relative assessment of the perfusion at the monitoring site. PI is displayed on a 10 segment LED bar, ranging from <1% (very weak perfusion) to >5% (strong perfusion). The highest LED will remain lit continuously to allow a PI level to be viewed. The Perfusion Index is the ratio of the AC (pulsatile) to DC (non-pulsatile) components of the IR (Infrared) signal where the AC and DC components correspond to the pulsatile and non-pulsatile amounts of blood, respectively.
⑭ Low SIQ	Flashes to indicate low SpO ₂ Signal IQ. Refer to Section 4, Low Signal IQ.
⑮ Alarm Silenced Indicator	Flashes to indicate the alarm is temporarily silenced.
⑯ Alarm Silence Button	Push once to temporarily silence the alarm for 120 seconds. Push a second time to return the unit to standard alarm monitoring.
⑰ Up button Down button	During saturation monitoring, use these buttons to adjust the volume of the pulse beep tone. Within the menu/setup system, these buttons are used to select values within each menu option.
⑱ Speaker	Provides audible indication of alarm conditions, pulse tone and feedback for key-presses. Ensure the speaker is not covered or the unit is placed face-down on bedding or other sound absorbing surface.

Rad-57cm Front Panel Controls




Rad-57cm

CONTROL / INDICATOR	DESCRIPTION
① Patient Cable Connector	Connects to a Rainbow sensor, Red sensor or a Red Patient Cable with an LNOP, LNOPv or LNCS sensor.
② Saturation (%SpO ₂) and Carboxyhemoglobin (%SpCO) Displays	The functional arterial hemoglobin oxygen saturation is displayed in units of SpO ₂ . When searching for a saturation and pulse, it will flash dashed lines. See <i>SpCO/SpMet Button</i> description, below.
③ Low PI Indicator	The Perfusion Index (PI) LED illuminates a constant light when arterial pulsations are very low (weak perfusions).
④ %SpCO Bar	Illuminates when SpCO capable sensor is attached. Bar will flash for SpCO alarm conditions. Continuously indicates the concentration of carboxyhemoglobin in 3% increments.
⑤ Pulse Rate Display and Methemoglobin (%SpMet) Displays	The pulse rate in beats per minute (bpm). When searching for a saturation and pulse, it will flash dashed lines. See <i>SpCO/SpMet</i> button description below

CONTROL / INDICATOR	DESCRIPTION
⑥ SpCO/SpMet Button	Pressing this button will display the numeric SpCO value in place of the SpO ₂ numeric value and the SpMet value in place of the pulse rate. Pressing the Mode/Enter or Next button will return to the SpO ₂ and pulse rate numeric value.
⑦ Mode / Enter Button	Used to enter the setup menus and to select/activate certain entries within the menu/setup system.
⑧ Next Button	Used within the menu/setup system to move through setup options. Not active during normal patient monitoring.
⑨ Power On / Off	Press to turn the unit on. Press-and-hold for 2 seconds to turn the unit off.
⑩ Battery Level Indicator	Four LED's indicate the status of the battery. When the final indicator begins flashing, replace the batteries.
⑪ Visual Alarm Indicator	Illuminates when any alarm condition exists. This indicator may not be turned-off or otherwise over-ridden.
⑫ %SpO ₂ / %SpCO Indicator	Indicator above label will illuminate to provide an additional visual indication of the value currently being displayed.
⑬ %SpMet Bar	Bar will flash for SpMet alarm conditions. Continuously indicates the concentration of Methemoglobin in increments: .5, 1-5, 7.5, 10, 15 >20%.
⑭ %SpMet/BPM	Indicator above label will illuminate to provide an additional visual indication of the value currently being displayed.
⑮ Low SIQ	Flashes to indicate low SpO ₂ Signal IQ. Refer to Section 4, <i>Low Signal IQ</i> .
⑯ Alarm Silenced Indicator	Flashes to indicate the alarm is temporarily silenced.
⑰ Alarm Silence Button	Push once to temporarily silence the alarm for 120 seconds. Push a second time to return the unit to standard alarm monitoring.
⑱ Up button Down button	During saturation monitoring, use these buttons to adjust the volume of the pulse beep tone. Within the menu/setup system, these buttons are used to select values within each menu option.
⑲ Speaker	Provides audible indication of alarm conditions, pulse tone and feedback for key-presses. Ensure the speaker is not covered or the unit is placed face-down on bedding or other sound absorbing surface.

Rad-57 Rear Panel (all models)



CONTROL / INDICATOR	DESCRIPTION
①	Serial Number Label Located on outside of case
②	Agency Approvals Label Located on outside of case
③	Battery Cover Located on back of unit
④	Battery Cover Release Press down and slide the battery cover off the bottom of the oximeter

SYMBOLS

SYMBOL	DESCRIPTION
	Caution, consult accompanying documents
	Type BF applied part complying with IEC 60601-1
	WEEE Compliant

Introduction

Before the Rad-57 Pulse CO-Oximeter can be used in a clinical setting, it needs to be inspected, properly setup and the batteries need to be installed.

Unpacking and inspection

Remove the instrument from the shipping carton and examine for signs of shipping damage. Check all materials against the packing list. Save all packing materials, invoice and bill of lading. These may be required to process a claim with the carrier.

If anything is missing or damaged, contact the Technical Service Department. The contact address and phone numbers are listed in Section 9, *Service and Repair*.

Preparation for monitoring

The following sections of the manual describe the preparation, set-up and initial installation of the Rad-57 Pulse CO-Oximeter.

POWER REQUIREMENTS

The Rad-57 is powered by 4 "AA" alkaline batteries. Do not use any other type of batteries or power source to run the device. The battery compartment is accessed from the back of the device. To install the batteries first remove the battery cover by depressing the small rectangular button at the bottom of the cover, and sliding the cover down off the bottom of the device. Install the batteries in the directions indicated by the battery icons inside the battery compartment. Replace the battery cover by sliding it back up from the bottom of the device until the rectangular locking button snaps back into position.

WARNING: USE ONLY ALKALINE BATTERIES. USE OF NON ALKALINE BATTERIES MAY AFFECT THE ACCURACY OF THE BATTERY LEVEL METER.

WARNING: USE OF BATTERIES WITH A CELL VOLTAGE OF MORE THAN 1.5V COULD CAUSE DAMAGE TO THE RAD-57.

Battery charge level is indicated by four LED indicators at the bottom of the front panel. All four indicators will be lit when the batteries are full, with fewer indicators being lit as the batteries lose their charge. When less than ten (10) percent battery life remains, the final battery indicator will begin to flash and an audible alarm will sound.

Monitor Setup

INITIAL SETUP

1. Inspect the oximeter case for damage.
2. Install 4 (four) new AA alkaline batteries.
3. Turn the unit on, the LEDs will scroll in the display window as the sensor calibrates, verify all indicators illuminate and speaker sounds a brief tone.
4. Configure the unit for your regional power line frequency (50 or 60 hz) if needed. Default is 60 hz (standard for the United States). See Section 4, *Special Menu, Special menu - Line Frequency Configuration*.

CAUTION: THE UNIT MUST BE CONFIGURED TO MATCH YOUR LOCAL POWER LINE FREQUENCY TO ALLOW FOR THE CANCELLATION OF NOISE INTRODUCED BY FLUORESCENT LIGHTS AND OTHER SOURCES.

No other setup is required. Refer to Section 4, *General Setup and Use* for additional steps to verify proper functioning of the unit.

Introduction

To operate the Rad-57 Pulse CO-Oximeter effectively, the operator must:

- Know how the oximeter derives its readings (see Section 1, *Pulse CO-Oximetry*)
- Be familiar with its controls and operation.
- Understand its status and alarm messages (see Section 5, *Alarm Identification, System Messages* and Section 6, *Troubleshooting*).

Basic Operation

GENERAL SETUP AND USE

1. Inspect the oximeter case for damage.
2. Ensure that the batteries are correctly installed.
3. Connect a Rainbow Sensor or a Red Patient Cable with an LNOP, LNOPv or LNCS sensor to the Patient Cable connector of the oximeter. Make sure it is a secure connection and the cable is not twisted, sliced or frayed. See Section 5, *messages*, to view messages that may be displayed pertaining to sensors and cables.
4. Select a sensor that is compatible with the oximeter before connecting it to the patient cable. See Section 8, *Sensors and Patient Cables*. If using a single patient adhesive or disposable sensor, check that the emitter (red light) and the photodetector are properly aligned. If using a reusable sensor, make sure it opens and closes smoothly. Remove any substances that may interfere with the transmission of light between the sensor's light source and photodetector.
5. Attach the sensor to the patient. Refer to the Directions for Use of the sensor.
6. Properly align and insert the male-connector end of the sensor into the female-connector end of unit (or patient cable). Make sure it is a secure connection.
7. Press the Power button to turn the oximeter on.
8. Verify all front-panel indicators momentarily illuminate and a one-second tone is heard.
9. Verify the front panel display is free of alarm and system failure messages (see Section 5, *Alarms and Messages*) and the battery indicator shows sufficient charge (see Section 4, *Battery Level Indicator*)
10. On the display, for each corresponding parameter and applicable unit, verify the readings for SpO₂, SpCO, SpMet and pulse rate (See Model Summary for applicable unit).

NOTE: " - - " will flash on the numeric display until the SpO₂, SpCO, SpMet and pulse rate readings have stabilized (less than 20 seconds for SpO₂ and up to 25 seconds for SpCO and SpMet).
11. Verify that the patient alarms are functional by setting the high and low SpO₂ and pulse rate alarm limits beyond the patient readings.
 - An alarm tone sounds.
 - The violated alarm parameter flashes.

12. Verify the sensor alarms are functional by removing the sensor from the sensor site.
 - "SEn OFF" message appears on the display.
 - The alarm tone sounds.
 - The Visual Alarm Indicator flashes.
 - Disconnect the sensor from the patient cable or oximeter.
 - Confirm that "NO SEn" message appears on the display.

Note: "NO SEn" and "SEn OFF" will only generate an alarm if the Rad-57 was actively monitoring a patient when the sensor was disconnected.
13. Verify parameter violation alarm silence operation.
 - Create an alarm condition by lowering the SpO₂ or pulse rate high alarm limits beyond the patient readings.
 - Press the Alarm Silence button.
 - The alarm tone ceases for 120 seconds.
14. To begin patient monitoring:
 - Adjust the alarm limits.
 - Adjust the alarm volume.
 - Adjust the pulse beep volume.
15. Verify the sensor is applied correctly and that the measured data is appropriate, see Section 4, *Successful SpO₂ Monitoring*.
16. Monitor the patient.
17. After monitoring is complete, remove the sensor from the patient and store or dispose of the sensor according to governing rules. See the Directions for Use of the sensor.
18. Press and hold the Power On/Off button for 2 seconds to turn the oximeter off.

NOTE: Turn the oximeter off between patients. This allows the unit to re-calibrate in order to interpret new physiological data and to conserve battery life.

DEFAULT SETTINGS

The Rad-57 stores two types of default values: those that the device automatically reverts to after a power cycle, and those that can be changed by the user and will be remembered after a power cycle.

The following table outlines the default values that the Rad-57 reverts to after a power cycle:

OPTION	DEFAULT SETTING
SpO ₂ high alarm limit	Set to Off
SpO ₂ low alarm limit	Set to 90%
Pulse rate high alarm limit	Set to 140 BPM
Pulse rate low alarm limit	Set to 50 BPM
SpCO* high alarm limit	Set to pre-power down setting (Factory Default is 10%)
SpCO* low alarm limit	Set to pre-power down setting (Factory Default is off)
SpMet* high alarm limit	Set to pre-power down setting (Factory Default is 5%)
SpMet* low alarm limit	Set to pre-power down setting (Factory Default is off)
Averaging Time	Set to pre-power down setting
FastSat	Set to pre-power down setting
Sensitivity	Set to pre-power down setting
Display brightness	Set to pre-power down setting
Pulse tone volume	Set to pre-power down setting
Alarm Volume	Set to pre-power down setting
Line Frequency	Set to pre-power down setting
Trend Active	Set to pre-power down setting
Alarm Silence	Set to all alarms active

*For each corresponding parameter, refer to the applicable unit.

Successful SpO₂ Monitoring

The following general points will aid in ensuring oximetry monitoring success.

- Place the sensor on a site that is not too thick, has sufficient perfusion and provides proper alignment of the LED's and photodetector.
- Place the sensor on a site that has unrestricted blood flow.
- Do not secure a sensor with tape.
- Do not select a site near potential electrical interference (electrosurgical unit, for example).
- Read the sensor Directions for Use for proper sensor application.

NUMERIC DISPLAY - SpO₂

Stability of the SpO₂ readings may be a good indicator of signal validity. Although stability is a relative term, experience will provide a good feeling for changes that are artifactual or physiological and the speed, timing, and behavior of each. The stability of the readings over time is affected by the averaging mode being used. The longer the averaging time, the more stable the readings tend to become. This is due to a dampened response as the signal is averaged over a longer period of time than during shorter averaging times. However, longer averaging times delay the response of the oximeter and reduce the measured variations of SpO₂ and PR.

MASIMO SENSORS

Before use, carefully read the LNOP, LNOPv, LNCS and Rainbow series sensor Directions for Use.

Use only Masimo oximetry sensors for SpO₂ measurements.

Tissue damage can be caused by incorrect application or use of an LNOP, LNOPv, LNCS or Rainbow sensor, for example by wrapping the sensor too tightly. Inspect the sensor site as directed in the sensor Directions for Use to ensure skin integrity and correct positioning and adhesion of the sensor.

CAUTIONS

- DO NOT USE DAMAGED LNOP, LNOPv, LNCS OR RAINBOW SENSORS. DO NOT USE AN LNOP, LNOPv, LNCS OR RAINBOW SENSOR WITH EXPOSED OPTICAL OR ELECTRICAL COMPONENTS. DO NOT IMMERSER THE SENSOR IN WATER, SOLVENTS, OR CLEANING SOLUTIONS (THE SENSORS AND CONNECTORS ARE NOT WATERPROOF). DO NOT STERILIZE BY IRRADIATION, STEAM, AUTOCLAVE OR ETHYLENE OXIDE (UNLESS OTHERWISE INDICATED ON THE SENSOR DIRECTIONS FOR USE). SEE THE CLEANING INSTRUCTIONS IN THE DIRECTIONS FOR USE FOR REUSABLE MASIMO LNOP, LNOPv, LNCS OR RAINBOW SENSORS.
- DO NOT USE DAMAGED PATIENT CABLES. DO NOT IMMERSER THE PATIENT CABLES IN WATER, SOLVENTS, OR CLEANING SOLUTIONS (THE PATIENT CABLE CONNECTORS ARE NOT WATERPROOF). DO NOT STERILIZE BY IRRADIATION, STEAM, AUTOCLAVE OR ETHYLENE OXIDE.

NUMERIC DISPLAY - PULSE RATE

The Pulse Rate displayed on the Rad-57 may differ slightly from the heart rate displayed on ECG monitors due to differences in averaging times. There may also be a discrepancy between cardiac electrical activity and peripheral arterial pulsation. Significant differences may indicate a problem with the Signal IQ due to physiological changes in the patient or one of the instruments or application of the sensor or patient cable. The pulsations from intra-aortic balloon support can be additive to the pulse rate displayed on the Pulse CO-Oximeter.

NUMERIC DISPLAY - SpCO (for applicable models)

A stable SpCO reading is associated with correct sensor placement, small physiological changes during the measurement and acceptable levels of arterial perfusion at the measurement site. Physiological changes at the measurement site are mainly caused by fluctuations in the oxygen saturation, blood concentration and perfusion.

Inaccurate measurements may be caused by:

- Significant levels of methemoglobin.
- Intravascular dyes such as indocyanine green or methylene blue.
- Abnormal hemoglobin levels.
- Abnormally low arterial perfusion.

NUMERIC DISPLAY - SpMet (for applicable models)

A stable SpMet reading is associated with correct sensor placement, small physiological changes during the measurement and acceptable levels of arterial perfusion at the measurement site. Physiological changes at the measurement site are mainly caused by fluctuations in the oxygen saturation, blood concentration and perfusion.

Inaccurate measurements may be caused by:

- Intravascular dyes such as indocyanine green or methylene blue.
- Abnormal arterial perfusion

LOW SIGNAL IQ (LOW SIQ)

The Rad-57 display provides a visual indicator Signal IQ and an alert when the displayed SpO₂ values are not based on adequate Signal IQ. The Signal IQ indicator is displayed on the Rad-57 as "Low SIQ".

The Low SIQ indicator flashes when the SpO₂ measurement may be compromised. When the Low SIQ indicator is flashing, proceed with caution and do the following:

- Assess the patient.
- Check the sensor and ensure proper sensor application. The sensor must be well secured to the site for the Rad-57 to maintain accurate readings. Also, misalignment of the sensor's emitter and detector can result in smaller signals.
- Determine if an extreme change in the patient's physiology and blood flow at the

monitoring site occurred, (e.g. an inflated blood pressure cuff, a squeezing motion, sampling of an arterial blood specimen from the hand containing the pulse oximetry sensor, severe hypotension, peripheral vasoconstriction in response to hypothermia, medications, or a spell of Raynaud's syndrome.)

- With neonates or infants, check that the peripheral blood flow to the sensor site has not been interrupted. For example, as may occur while lifting or crossing their legs during a diaper change.

After performing the above, if the "Low SIC" indication occurs frequently or continuously, obtaining an arterial blood specimen for CO-Oximetry analysis may be considered to verify the oxygen saturation value.

LOW PERFUSION

Low Perfusion is indicated when the arterial pulsations are very low (weak perfusion). Low Perfusion is shown when the Low PI LED indicator light is illuminated constantly.

CAUTION: IF THE LOW PERFUSION INDICATION IS FREQUENTLY DISPLAYED, FIND A BETTER PERFUSED MONITORING SITE. IN THE INTERIM, ASSESS THE PATIENT AND, IF INDICATED, VERIFY OXYGENATION STATUS THROUGH OTHER MEANS.

ACTIONS TO BE TAKEN

If the SpO₂ readings show significant differences, do the following:

- Make sure the emitter and photodetector are aligned directly opposite each other.
- Select a site where the distance between the emitter and photodetector is minimized.
- Wipe the sensor site with a 70% isopropyl alcohol pad or rubefacient cream (10-30% methyl salicylate and 2-10% menthol) and allow to dry for 20-30 seconds. Strong vasodilator creams, such as nitroglycerin paste, are not recommended.
- If possible, remove electrical noise sources such as electrosurgical units or other electrical/electronic equipment.
- If artificial nails or excessive fingernail polish are present, select another site or remove the polish/artificial nails.
- If possible, ensure that the sensor is placed in a location with low ambient light.

CAUTION: IF ANY MEASUREMENT SEEMS QUESTIONABLE, FIRST CHECK THE PATIENT'S VITAL SIGNS BY ALTERNATE MEANS AND THEN CHECK THE PULSE CO-OXIMETER FOR PROPER FUNCTIONING.

SENSITIVITY

Three sensitivity levels enables a clinician to tailor the response of the Rad-57 to the needs of the particular patient situation. They are as follows:

- Normal Sensitivity – This is the recommended mode for typical monitoring purposes. It is advisable for care areas where patients are observed frequently, such as ICUs.

■ Adaptive Probe Off Detection (APOD™) – This is the recommended monitoring mode where there is a high probability of the sensor becoming detached. It is also the suggested mode for care areas where patients are not visually monitored continuously. This mode delivers enhanced protection against erroneous pulse rate and arterial oxygen saturation readings when a sensor becomes inadvertently detached from a patient due to excessive movement.

■ Maximum Sensitivity (MAX) - This mode is recommended for patients with low perfusion or when the low perfusion message is displayed on the screen in APOD or normal sensitivity mode. This mode is not recommended for care areas where patients are not monitored visually, such as general wards. It is designed to interpret and display data at the measuring site when the signal may be weak due to decreased perfusion. When a sensor becomes detached from a patient, it will have compromised protection against erroneous pulse rate and arterial saturation readings

CAUTION: When using the Maximum Sensitivity setting, the performance of the **SENSOR OFF** detection may be compromised. If the unit is in this setting and the sensor becomes dislodged from the patient, the potential for false readings may occur due to environmental 'noise' such as light, vibration and excessive air movement.

BATTERY LEVEL INDICATOR

Four LED indicators provide information on the remaining battery capacity. The operator should monitor these indicators periodically to determine remaining battery life and if the batteries should be replaced. Battery capacity is indicated in the following chart.

INDICATION	BATTERY CAPACITY
4 LEDs	100% to 75%
3 LEDs	75% to 50%
2 LEDs	50% to 25%
1 LED	25% to 10%
1 FLASHING LED WITH AUDIBLE ALARM	10% to 0%

LOW BATTERY AUDIBLE ALARM

If a low battery condition occurs during patient monitoring, a low priority alarm will sound, and can be silenced for 120 seconds by pressing the Alarm Silence Button.

If a low battery condition occurs while not monitoring a patient, pressing the Alarm Silence Button will suspend the the alarm until the power is cycled or patient monitoring begins.

If a low battery condition occurs, immediately discontinue patient monitoring and replace the batteries.

NOTE: Remove batteries when storing unit for prolonged periods to maintain battery life.

WARNING: FAILURE TO REPLACE BATTERIES PROMPTLY AFTER A LOW BATTERY ALARM MAY RESULT IN THE CO-OXIMETER SHUTTING DOWN AND LEAVING THE PATIENT IN AN UNMONITORED CONDITION.

WARNING: USE ONLY ALKALINE BATTERIES. USE OF NON ALKALINE BATTERIES MAY AFFECT THE ACCURACY OF THE BATTERY LEVEL METER.

WARNING: USE OF BATTERIES WITH A CELL VOLTAGE OF MORE THAN 1.5V COULD CAUSE DAMAGE TO THE RAD-57.

WARNING: EFFECTIVE BATTERY LIFE WILL BE REDUCED WHEN OPERATING THE INSTRUMENT BELOW 5°F (-15° C) DUE TO ALKALINE BATTERY TECHNOLOGY.

Normal Patient Monitoring








During normal operation, the Rad-57 Display shows oxygen saturation (as % SpO₂) on the upper number and Pulse Rate (in beats per minute) on the lower number.

On each model, a 10-LED bar displays the following parameters:

- Carboxyhemoglobin saturation (%SpCO) and Perfusion Index (PI) on the Rad-57c.
- Methemoglobin saturation (%SpMet) and Perfusion Index (PI) on the Rad-57m.
- Carboxyhemoglobin saturation (%SpCO) and Methemoglobin saturation (%SpMet) on the Rad-57cm. Low Perfusion Index (PI) is displayed with an LED on this model.

The following sections describe the function of the Rad-57 front panel controls during normal patient monitoring.

RAD-57 FRONT PANEL CONTROL OPERATION

BUTTON	FUNCTION
 SpCO Rad-57c	Rad-57c and Rad-57cm: Pressing the applicable button will display the numeric SpCO/SpMet value for 10 seconds in place of the SpO ₂ or pulse rate numeric value. Pressing the "Mode/Enter" or "Next" button during this 10-second period will return to the SpO ₂ numeric value.
 SpMet Rad-57m	Rad-57m: Pressing the applicable button will display the numeric SpMet value in place of the pulse rate numeric value. Pressing the button again will return to displaying the pulse rate.
 MODE ENTER	Enters the Rad-57 setup/menu system. See Section 4, <i>Operation</i> .
 NEXT	No function during normal patient monitoring.
	Push once to temporarily silence the alarm for 120 seconds. Push a second time to return the unit to standard alarm monitoring.
	During normal patient monitoring, the "Up" and "Down" Arrow keys control the Pulse Tone volume. At the lowest setting, the pulse tone is muted. A low-pitch tone indicates the highest or lowest setting has been reached. In the setup/menu system, the "Up" and "Down" Arrow keys select among the options for each setting.
	Power "on/off" button. Press this button to turn the unit on. Press-and-hold for 2 seconds to turn the unit off.

Setup Menu



This section gives an overview of the Rad-57 menu selections available. To navigate through the menus, use the *Mode/Enter*, *Next*, *Up* and *Down* keys located on the front panel of the Pulse CO-Oximeter, below the LED display. The following sub sections describe each menu item in more detail. The Pulse CO-Oximeter has options that allow user configuration to suit specific needs.

MENU NAVIGATION

The Rad-57 set-up and configuration options are accessed through the menu system. The *Mode/Enter* key is used to enter the menu system and to move through the different menu levels. Within each level of the system, the *Next* key is used to move from one option to the next. The *Up* and *Down* arrow keys are used to select values within each option. The parameter is set/selected when either the *Mode/Enter* or *Next* keys are pressed.

SETUP MENU LEVEL 1 – ALARM VOLUME









Push the *Mode/Enter* button to enter menu level 1.

SETTING	
 3X	Alarm Volume
	Alarm on/off

Use *Up* or *Down* Arrow Keys to adjust parameter to desired setting.

SETUP MENU LEVEL 2 – ALARM LIMITS

Push the *Mode/Enter* button again to enter menu level 2.




SETTING	
 2X	SpO ₂ Low Alarm Limit
	SpO ₂ High Alarm Limit
	Pulse Rate Low Alarm Limit
	Pulse Rate High Alarm Limit
	SpCO* Low Alarm Limit
	SpCO* High Alarm Limit
	SpMet* Low Alarm Limit
	SpMet* High Alarm Limit

Use *Up* or *Down* Arrow Keys to adjust parameter to desired setting.

*See Model Summary for applicable unit.

SETUP MENU LEVEL 3 – AVERAGING AND SENSITIVITY

Push the *Mode/Enter* button again to enter menu level 3.

SETTING	
 3X	Sensitivity HI = Maximum Nor = Normal APO = APOD
	Averaging. The signal averaging time of this device can be set to: 2*, 4*, 8, 10, 12, 14 or 16 seconds
	FastSat* On, Off

Use *Up* or *Down* Arrow Keys to adjust parameter to desired setting.

Note: These changes affect SpO₂ monitoring only.














*Select "Yes" to activate the FastSat algorithm. The FastSat averaging time is dependant on the input signal. FastSat is automatically enabled in 2 and 4 second averaging.

SETUP MENU LEVEL 4 - TREND SETTINGS

Push the *Mode/Enter* button again to enter menu level 4.

To enable trending of patient data, the trend feature must be enabled (set to ON), and the current date and time must be set. See Section 4, *Trend setup and use*.

The current date and time can only be set if the Trend is set to "ON". The date and time menu selections are not available if Trend is set to "OFF".

SETTING	
 4X	Trend ON / OFF
	Sequence for Rad-57c and Rad-57cm
	Set Month
	Set Day
	Set Year
	Set Hour
	Set Minute
	Sequence for Rad-57m
	Set Year
	Set Month
	Set Day
	Set Hour
	Set Minute

Use *Up* key to turn trend ON. Use *Down* key to turn trend OFF.

Use *Up* or *Down* Arrow Keys to adjust parameter to desired setting.



A valid date must be entered. If an invalid date is entered (i.e. February 31), the trend will not turn on and "In d off" will be displayed.

Note: The date and time must be set before trending will be enabled. The Rad-57 will automatically 'time out' of the setup menu after 10 seconds with no key presses. If the Rad-57 should time-out of the Trend Settings menu, the trend will not be enabled.

Note: Enabling trend (setting Trend to "ON") will erase all trend information in the Rad-57

SETUP MENU LEVEL 5 - LED BRIGHTNESS AND FACTORY DEFAULTS

Push the Mode/Enter button again to enter menu level 5.

SETTING	
	LED Display Brightness (4 levels) Note: All LED indicators are illuminated while adjusting this setting.
	Restore Factory Defaults
Use Up or Down Arrow Keys to adjust parameter to desired setting.	

Pressing  a sixth time returns the Rad-57 to patient monitoring display in the Saturation/Pulse Rate Mode.

MENU SELECTION

SpO₂ and pulse rate will be continuously displayed once the sensor is properly placed and calibration is complete. To display continuous SpCO and SpMet monitoring (see Model Summary for applicable unit), press . To return to continuous SpO₂ monitoring and pulse rate view, press  again.

Note: When in continuous SpCO/SpMet monitoring mode and a Multi-Parameter Alarm Condition occurs, the instrument will automatically exit continuous SpCO/SpMet monitoring mode and revert to SpO₂ monitoring mode.

POWER OFF

OFF – Press and hold the On/Off button for two seconds.





Special Menu



This section gives an overview of the Rad-57 special menu selection available. To navigate through the menu, use the Power, Next, Up and Down keys located on the front panel of the oximeter. The oximeter has options that allow user configuration to suit specific needs.

The following sub-section describes this menu item in more detail.

SPECIAL MENU – LINE FREQUENCY CONFIGURATION

NOTE: The unit must be configured to match your local power line frequency (50 or 60 hz) to operate properly. The unit default is set to 60 hz (standard for the United States).

1. Turn the Rad-57 off.
2. Hold the Down  arrow key while turning the Rad-57 back on.
3. Push the Next  button 5 times. "LF" will be displayed in the top LED display window and the active line frequency will be displayed in the lower LED display window.
4. Push the Up  arrow to set the line frequency to 60 Hz and the Down  arrow to set it to 50 Hz.
5. Turn the unit off.

SETTING	
Hold  + press 	Enter Line Frequency Menu
LF	Set Line Frequency
Use Up or Down Arrow Keys to adjust parameter to desired setting. Note: The parameter is set/selected when the unit is turned off.	

Trend Setup and Use

INTRODUCTION

The Rad-57 can store 72 hours of SpO₂, Pulse Rate, SpCO*, SpMet*, and Perfusion Index trend data, captured at 2 second intervals. This trend data can then be transferred to a PC for evaluation.

Trend data is stored in non-volatile memory, so it is not erased when the unit is shut off or when the batteries are replaced.

A Pronto serial cable is required to connect the sensor connector of the Rad-57 to the PC and TrendCom software is required to download the trend data on to the PC. Refer to Section 10 - Accessories and contact Masimo to acquire these optional items. Patient monitoring is not possible while trend memory is being transferred to a PC.

A trend data download is initiated using the TrendCom utility which downloads the trend data and saves it to a space-delimited ASCII text (.out) file.

*See Model Summary for applicable unit.

TRENDCOM UTILITY INSTALLATION

Copy the TrendCom utility from the CD onto a PC running MS-Windows.

TRENDCOM UTILITY OPERATION

1. Disconnect patient cable from the Rad-57.
2. Connect the mini-D end of the Rad-57 PRONTO serial cable to the Rad-57 patient cable connector (See Section 2, *Rad-57 front panel controls*) and connect the DB-9 end to a COM port on the PC.
3. Turn the Rad-57 on
4. Start the TrendCom Utility
5. Select the appropriate COM port number, if necessary.
6. Push the **RETRIEVE TREND** button on the TrendCom utility.
Select the desired location and assign a filename for the trend file. Press **Save**.
7. The Rad-57 will display "dat out" while trend data is being transferred.
A progress bar will advance to indicate the status of the download.
Larger trend files will take longer to download. Transfer time is approximately 20 seconds per hour of trend data.
Note: *During download of trend information, all normal Rad-57 functions are unavailable and the keypad is locked, except for the power button.*
8. When trend data transfer is complete, close TrendCom and disconnect the Rad-57 from the Rad-57 PRONTO serial cable.
9. Turn the Rad-57 off to exit the trend download mode.
Note: *USB to serial port adapters are not supported for trend transfer.*
Note: *Enabling trend (setting Trend to "ON") will erase all trend information in the Rad-57.*

ERASING TREND MEMORY

To erase (clear) the trend memory, turn the trend off and back on again. Enabling trend (setting Trend to "ON") will erase all trend data.

Note: *Turning trend off will not erase trend memory. You may turn trending off and still retrieve the trend data using TrendCom.*

Turning the Rad-57 off or replacing the batteries will not erase the trend data.

Turn trending off before storing the unit for any length of time.

TREND DATA FORMAT

After a successful download of the trend data, a .out file will be created containing the trend-dump information in ASCII delimited format. The format is defined in the following table.

PARAMETER	SPECIFICATION
Date	MM/DD/YY
Time	HH:MM:SS
SpO ₂ *	001 to 100, or "---" meaning parameter not available
SpCO*	00.00 to 100, or "---" meaning parameter not available
SpMet*	00.00 to 100, or "---" meaning parameter not available
Pulse Rate	001 to 240, or "---" meaning parameter not available
Perfusion Index	00.00 to 20.00
The exceptions are displayed as a 3 digit, ASCII encoded, hexadecimal value. The binary bits of the hexadecimal value are encoded as follows:	
	000 = Normal operation; no exceptions
	001 = No Sensor
	002 = Defective Sensor
	004 = Low Perfusion
	008 = Pulse Search
	010 = Interference
	020 = Sensor Off
	040 = Ambient Light
	080 = Unrecognized Sensor
	100 = reserved
	200 = reserved
	400 = Low Signal IQ
	800 = Masimo SET. This flag means the algorithm is running in full SET mode. It requires a Masimo sensor and needs to acquire some clean data for this flag to be set
Exception Messages	

*See Model Summary for applicable unit.

SAMPLE TREND OUTPUT

```

07/21/04 09:56:08 SpO2=000 SpCO=000 PR=000 PI=00.00 EXC=820:OffPat,SET
07/21/04 09:56:10 SpO2=000 SpCO=000 PR=000 PI=00.00 EXC=828:Search,OffPat,SET
07/21/04 09:56:12 SpO2=097 SpCO=001 PR=069 PI=04.69 EXC=800:SET
07/21/04 09:56:14 SpO2=096 SpCO=001 PR=074 PI=02.28 EXC=C00:LowSigIO,SET
07/21/04 09:56:16 SpO2=098 SpCO=001 PR=078 PI=03.64 EXC=800:SET
07/21/04 09:56:18 SpO2=000 SpCO=000 PR=000 PI=00.00 EXC=800:SET
07/21/04 09:56:20 SpO2=000 SpCO=000 PR=000 PI=00.00 EXC=820:OffPat,SET
07/21/04 09:56:22 SpO2=096 SpCO=001 PR=078 PI=02.68 EXC=800:SET

```

Alarm Indication

An alarm condition is indicated by:

- Audible alarm tone
 - Visual Alarm Indicator
 - Out-of-limit parameter will flash
- "SEn OFF" and "nO SEN" will only generate an alarm condition after a pulse has been found.

Alarm Limits

CAUTION: TO ENSURE THAT ALARM LIMITS ARE APPROPRIATE FOR THE PATIENT BEING MONITORED, CHECK THE LIMITS EACH TIME THE PULSE CO-OXIMETER IS USED.

An audible alarm and a flashing alarm status indicator will occur when an alarm limit is met or exceeded for greater than five seconds. It is best that the operator be within a minimum of 10 feet from the unit. Directions for alarm suspension are indicated below. When a sensor is not connected to a patient, or when a sensor is not connected to its cable, the display will read SEn OFF or NO SEN. An audible alarm will accompany the display unless the oximeter has been set to Alarm Silence Mode.

SETTING	RANGE
SpO ₂ High Limit	The SpO ₂ high alarm limit can be set anywhere between 2% and 100%, with a 1% step size. In the "----" (off) setting, the SpO ₂ High Limit alarm is disabled.
SpO ₂ Low Limit	The SpO ₂ low alarm limit can be set anywhere between 1% and 100%, with a 1% step size. Note: The low alarm limit must always be below the high alarm limit. If the high alarm limit is set below the low alarm limit, the low alarm limit will automatically adjust to the next setting below the newly entered high alarm limit.
Pulse Rate High Limit (BPM)	The pulse rate high alarm limit can be set anywhere between 30 BPM and 240 BPM, with a 5 BPM step size.
Pulse Rate Low Limit (BPM)	The pulse rate low alarm limit can be set anywhere between 25 BPM and 235 BPM, with a 5 BPM step size. Note: The low alarm limit must always be below the high alarm limit. If the high alarm limit is set below the low alarm limit, the low alarm limit will automatically adjust to the next setting below the newly entered high alarm limit.
SpCO* High Limit	The SpCO high alarm limit can be set anywhere between 5% and 50%, with a 5% step size. In the "----" (off) setting, the SpCO High Limit Alarm is disabled. Factory default setting is 10%.
SpCO* Low Limit	The SpCO low alarm limit can be set anywhere between 5% and 45%, with a 5% step size. In the "----" (off) setting, the SpCO Low Limit Alarm is disabled. Factory default setting is "off". Note: The low alarm limit must always be below the high alarm limit. If the high alarm limit is set below the low alarm limit, the low alarm limit will automatically adjust to the next setting below the newly entered high alarm limit.

*See Model Summary for applicable unit.

ALARM LIMITS CONTINUED

SETTING	RANGE
SpMet* High Limit	The SpMet high alarm limit can be set anywhere between 1% and 70% with a 1% step size. In the "....." (off) setting, the SpMet High Alarm Limit Alarm is disabled. Factory default setting is 5%.
SpMet* Low Limit	The SpMet low alarm limit can be set anywhere between .5% and 10% with a .5% step size. In the "....." (off) setting, the SpMet Low Limit Alarm is disabled. Factory default setting is "off". <i>Note: The low alarm limit must always be below the high alarm limit. If the high alarm limit is set below the low alarm limit, the low alarm limit will automatically adjust to the next setting below the newly entered high alarm limit.</i>

*See Model Summary for applicable unit.

Note: If there is a loss of power for any length of time, the Alarm settings will be set back to the User set defaults. If the user has not utilized this option, then they will be set back to the factory defaults.

ALARM SILENCE

Audible alarms may be suspended, while visual alarms may not. The alarm suspension setting is controlled by the Alarm Silence Button. Pressing the Alarm Silence Button will suspend the alarm for 120 seconds.

Power-On – Alarms are active and Alarm Silenced Indicator is off.

Push Once – Alarm is silenced for 120 seconds and Alarm Silenced Indicator flashes.

Push Twice - Return to Audible Alarm Active.

ALARM SILENCED INDICATOR

The Alarm Silenced Indicator provides visual feedback regarding the audible alarm status. The audible alarms are muted when the indicator is flashing.

While monitoring a patient, acknowledging an alarm condition by pressing the Alarm Silence Button (one time) will silence the alarm tone for 120 seconds and the Alarm Silenced Indicator will flash. Pressing the Alarm Silence Button a second time (while the Alarm Silenced Indicator is still flashing) will activate alarms and alarm indicator is off.

While not monitoring a patient, acknowledging an alarm condition by pressing the Alarm Silence Button (one or more times) will permanently silence the alarm tone, and the Alarm Silenced Indicator will remain illuminated until the power is cycled or patient monitoring begins.

Should the alarm condition be created by low batteries, replace the batteries before monitoring begins or continues.

MESSAGES

The Rad-57 will indicate other data or system errors.

Message conditions for the Rad-57 follow:

DISPLAY	TYPE	SOLUTION
SpO ₂ NUMBER FLASHES	Saturation limit alarm	Assess /address patient condition. Re-set alarm limits if indicated.
PULSE RATE NUMBER FLASHES	Pulse Rate limit alarm	Assess /address patient condition. Re-set alarm limits if indicated.
NO SEN	No Sensor Connected	Connect sensor to cable.
SEN OFF	Sensor off patient	1. Reattach sensor to patient. 2. Verify proper sensor placement.
02 SEN	SpO ₂ sensor attached	If SpCO and/or SpMet parameters are desired, attach a Rainbow sensor to the unit.
LEDS FLASH HORIZONTAL BARS	Pulse Search	Wait for pulse detection. (This search should occur whenever a sensor is first applied to a patient). If necessary, shield the sensor from excessive ambient or strobing light.
CIRCULATING LEDES	Sensor is calibrating	
LOW SIG INDICATOR FLASHES	Low SpO ₂ Signal IQ	1. Rule out occlusion of blood flow. 2. Verify placement of sensor.
PERFUSION INDEX (PI) BAR* TURNS RED (Bottom two LEDs only) OR LOW PI LED IS FLASHING	Low Signal Strength	1. Rule out occlusion of blood flow. 2. Attempt to warm patient. 3. Move sensor to better perfused site. Note: <i>Maximo recommends using an adhesive sensor whenever low perfusion is expected or evident.</i>
% SPCO* BAR GRAPH	Alarm Condition	Continuously indicates the concentration of Carboxyhemoglobin in 3% - 5% increments. When an SpCO alarm is violated, the bargraph will flash the color of the SpCO level being violated along with an audible alarm. Note: <i>For a more accurate reading press the SpCO button.</i>

*See Model Summary for applicable unit.

MESSAGES CONTINUED

DISPLAY	TYPE	SOLUTION
% SpMet* BAR GRAPH	Alarm Condition	Continuously indicates the concentration of Methemoglobin in the following increments: .5, 1-5, 7.5, 10, 15 and >20%. When an SpMet alarm is violated, the bargraph will flash at the SpMet level being violated along with an audible alarm.
SINGLE BATTERY LEVEL INDICATOR FLASHES (WITH AUDI- BLE ALARM)	Battery level too low	Replace batteries immediately.
NO CBL	No Cable Connected	Connect appropriate cable to unit.
Err ##	System Fault	Return for service. There are several error codes. All error codes require return of the unit to an authorized service center for repair. See Section 9, <i>Service and Repair</i> for return procedure.
bAd CBL	Defective cable	Replace cable
CBL (Blinking)	Incompatible cable	Connect appropriate cable
bAd SEN	Defective sensor	Replace sensor
SEN (Blinking)	Unrecognized sensor	Connect appropriate cable
INT DET (Blinking)	Interference detected	Ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required.

*See Model Summary for applicable unit.

Troubleshooting

The following chart describes what to do if the Rad-57 system does not operate properly or fails.

PROBLEM	POSSIBLE CAUSE(S)	RECOMMENDATION
DIFFICULTY OR NO SpCO₂/SpMet* READING	<p>Interference from line-frequency induced noise.</p> <p>Inappropriate sensor</p> <p>Excessive motion</p> <p>Excessive ambient or strobing light</p>	<p>Minimize or eliminate interference from surgical or fluorescent lighting.</p> <p>Verify/set 50/60hz menu setting. Refer to Section 3, <i>Initial Setup</i> for details.</p> <p>Verify use of an SpCO₂/SpMet capable sensor. Red LNOP, LNCS or LNOPv sensors can not provide SpCO₂ or SpMet measurements.</p> <p>Minimize or eliminate motion at the monitoring site.</p> <p>Also, see Section 4, <i>Successful Monitoring</i> for additional information.</p> <p>Shield the sensor from excessive light.</p>
UNIT DOES NOT POWER ON	<p>Low battery</p>	<p>Check / replace batteries</p>
CONTINUOUS SPEAKER TONE	<p>Internal Failure</p>	<p>Unit requires service. Press the Alarm Silence button. If alarm continues to sound, power down unit and remove batteries.</p>
NO SPEAKER TONE	<p>Pulse tone set to "mute"</p> <p>Alarm Silence Enabled</p>	<p>Press Up Arrow.</p> <p>Inspect Alarm Silence Indicator. See Section 4, <i>Alarm Silence</i>. Press Alarm Silence Button until Alarm Silence Indicator is no longer flashing.</p>
BUTTONS DON'T WORK WHEN PRESSED	<p>Internal Failure</p>	<p>Return for service.</p>
LOW BATTERY ALARM SOUNDS. BATTERY INDICATOR SHOWS LOW BATTERY CAPACITY LESS THAN EXPECTED CAPACITY	<p>Effective battery life will be reduced when operating the instrument below 5 degrees Fahrenheit due to alkaline battery technology.</p>	<p>Remove the batteries and allow them to warm up to room temperature, re-install them and check the battery indicator level. If the battery capacity remains low, replace batteries.</p>

*See Model Summary for applicable unit.

Rad-57 Family Specifications

PERFORMANCE

measurement range	1-100%
Oxygen Saturation (%SpO ₂):	1 - 99%
Carboxyhemoglobin Saturation (%SpCO*)	1 - 99.9%
Methemoglobin Saturation (%SpMet*)	25-240 beats per minute
Pulse Rate (bpm):	0.02 - 20%
Perfusion Index:	<1 second delay
Response time:	

ACCURACY

Saturation	70% to 100%
-------------------	-------------

*No Motion*¹

Adults, Pediatrics	±2 digits
Neonate	±3 digits

*Motion*²

Adults, Pediatrics	±3 digits
Neonate	±3 digits

*Low Perfusion*³

Adults, Pediatrics	±2 digits
Neonate	±3 digits

Carboxyhemoglobin Saturation Accuracy (%SpCO*)⁷

Adults, Pediatrics, Neonate	1% - 40% ±3 digits
-----------------------------	--------------------

Methemoglobin Saturation Accuracy (%SpMet)⁸

Adults, Pediatrics, Neonate	1% - 15% ±1 digit
-----------------------------	-------------------

Pulse Rate (bpm) during no motion conditions¹

Adults, Pediatrics, Neonate	25 to 240 ±3 digits
-----------------------------	---------------------

Pulse Rate (bpm) during motion conditions²

Adults, Pediatrics, Neonate	25 to 240 ±5 digits
-----------------------------	---------------------

Resolution

Oxygen Saturation (%SpO ₂)	1%
Carboxyhemoglobin saturation (%SpCO*), digital display	1%
Carboxyhemoglobin saturation(%SpCO*), continuous bar display	5%
Methemoglobin saturation (%SpMet*), digital display	.1%
Methemoglobin saturation(%SpMet*), continuous bar display	.5, 1-5, 7.5, 10, 15 >20%
Pulse Rate (bpm) ⁹	1 bpm

*See Model Summary for applicable unit.

Interfering Substances

Carboxyhemoglobin and methemoglobin may erroneously increase oxygen saturation readings. The level of increase is approximately equal to the amount of carboxyhemoglobin and/or methemoglobin present. Dyes, or any substance containing dyes, that change usual arterial pigmentation may cause erroneous readings.

ELECTRICAL**Batteries**

Type:	4 "AA" Alkaline ⁵
Capacity:	over 8 hours ⁴
Isolation:	No external power or ground connection, internally powered only, DC current.

ENVIRONMENTAL

Operating Temperature:	0°F to 129°F (-18°C to 54°C)
Storage Temperature:	-40°F to 158°F (-40°C to +70°C) ⁵
Operating Humidity:	5% to 95%, non-condensing
Operating Altitude:	500 mbar to 1060 mbar pressure, -1000 ft to 18,000 ft (-304 m to 5,486 m)

PHYSICAL CHARACTERISTICS

Dimensions:	6.2" x 3.0" x 1.4" (15.8 cm x 7.6cm x 3.6 cm)
Weight:	13oz. (0.37 kg)

SpO₂ Modes

Averaging mode:	2, 4, 8, 10, 12, 14 or 16 seconds ⁶
Sensitivity:	Normal, Maximum and APOD

Alarms*

Audible and visual alarms for high low saturation, pulse rate, carboxyhemoglobin and methemoglobin (SpO₂ range 1-100%, pulse rate range 25-240 bpm, SpCO 5%-50%, SpMet 1% - 70%)
Sensor condition, system failure and low battery alarms

High Priority:	571 Hz tone, 5 pulse burst, pulse spacing: 0.250s, 0.250s, 0.500s, 0.500s, 0.250s, repeat time: 10s
Low Priority:	500Hz tone, 1 pulse, repeat time: 5s
Volume:	75 dB (max.)

Display/Indicators*

Data display: %SpO₂, %SpCO, SpCO bar, pulse rate, %SpMet, SpMet bar, alarm status, alarm silenced status, Low Signal IQ, battery status.

Type:	LED
Display update rate	1 second

*See Model Summary for applicable unit.

Compliance

EMC Compliance:	EN60601-1-2, Class B
Equipment Classification:	IEC 60601-1
Type of Protection:	Internally powered (on battery power)
Degree of Protection-Patient Cable:	Type BF-Applied Part
Mode of Operation:	Continuous

- Masimo Rainbow SET technology with LNOP, LNOPV and LNCS sensors has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory CO-Oximeter and ECG monitor. 1% was added to account for the properties of fetal hemoglobin. This variation equals plus or minus one standard deviation which encompasses 68% of the population.
- Masimo Rainbow SET technology with LNOP, LNOPV and LNCS sensors has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory co-oximeter and ECG monitor. 1% was added to account for the properties of fetal hemoglobin. This variation equals plus or minus one standard deviation which encompasses 68% of the population.
- Masimo Rainbow SET technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and a % transmission of greater than 5% for saturations ranging from 70 to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.
- This represents approximate run time at lowest indicator brightness, using new, fully charged batteries.
- If alkaline batteries are to be stored for extended periods of time, it is recommended that they be stored between -0°C to +30°C, and at a relative humidity less than 85%. If stored for a prolonged period at environmental conditions beyond these limits, overall battery capacity may be diminished, and lifetime of the batteries may be shortened.
- With FastSat the averaging time is dependent on the input signal. For the 2 and 4 second settings the averaging time may range from 2.4 and 4-6 seconds, respectively.
- SpCO accuracy has been validated on healthy adult male and female volunteers with light to dark skin pigmentation in the range of 0% - 40% SpCO against a laboratory co-oximeter. This variation equals plus or minus one standard deviation which encompasses 68% of the population. The SpCO accuracy has not been validated under motion conditions.
- SpMet accuracy has been validated on healthy adult male and female volunteers with light to dark skin pigmentation in the range of 0% - 15% SpMet against a laboratory co-oximeter. The variation equals plus or minus one standard deviation which encompasses 68% of the population. The SpMet accuracy has not been validated under motion conditions.
- Masimo Rainbow SET technology with LNOP, LNOPV and LNCS sensors has been validated for pulse rate accuracy for the range of 25-240 bpm in bench top testing against a Biotek Index 2 simulator. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

Introduction

This section covers the use and cleaning of Masimo sensors and patient cables.

Before use of any sensor, carefully read the sensor's Directions for Use.

Use only Masimo oximetry sensors and cables for SpO₂, SpCO, SpMet measurements. Other oxygen transducers or sensors may cause improper Rad-57 Pulse CO-Oximeter performance.

Tissue damage can be caused by incorrect application or use of a sensor, for example by wrapping the sensor too tightly. Inspect the sensor site as directed in the sensor Directions for Use to ensure skin integrity and correct positioning and adhesion of the sensor.

CAUTIONS:

- DO NOT USE DAMAGED SENSORS OR PATIENT CABLES. DO NOT USE A SENSOR OR PATIENT CABLE WITH EXPOSED OPTICAL OR ELECTRICAL COMPONENTS.
- DO NOT IMMERGE THE SENSOR OR PATIENT CABLE IN WATER, SOLVENTS, OR CLEANING SOLUTIONS (THE SENSORS AND CONNECTORS ARE NOT WATERPROOF).
- DO NOT STERILIZE SENSORS OR PATIENT CABLES BY IRRADIATION, STEAM, AUTOCLAVE OR ETHYLENE OXIDE (UNLESS OTHERWISE INDICATED ON THE SENSOR DIRECTIONS FOR USE). SEE THE CLEANING INSTRUCTIONS IN THE DIRECTIONS FOR USE FOR REUSABLE MASIMO SENSORS.
- DO NOT REPROCESS ANY MASIMO SET PATIENT CABLES OR SINGLE USE SENSORS.
- ALL SENSORS AND CABLES ARE DESIGNED FOR USE WITH SPECIFIC MONITORS. VERIFY THE COMPATIBILITY OF THE MONITOR, CABLE AND SENSOR BEFORE USE, OTHERWISE PATIENT INJURY CAN RESULT.

SELECTING A SENSOR

When selecting a sensor, consider, the patient's weight, the adequacy of perfusion, the available sensor sites, and the duration of monitoring. For more information refer to the following table or contact your Sales Representative. Use only Masimo sensors and sensor cables. Select an appropriate sensor, apply it as directed, and observe all warnings and cautions presented in the Directions for Use accompanying the sensor.

High ambient and strobing light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor. To prevent interference from ambient or strobing light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.

SENSOR APPLICATION SITE

Unless indicated otherwise in the directions for use, reposition reusable sensors at least every 4 hours and for adhesive sensors inspect the site at least every 8 hours or sooner. If indicated by circulatory condition or skin integrity, reapply to a different monitoring site.

Masimo Rainbow® Sensors

Masimo Rainbow sensors must be used with the Rad-57 to enable measurement of Carboxyhemoglobin (SpCO) and Methemoglobin (SpMet). Rainbow sensors will only function with oximeter devices equipped with Masimo Rainbow SET Technology.

Rainbow sensors connect to the device directly or with a patient cable.

RAINBOW REUSABLE SENSORS

SpO₂ and pulse rate accuracy for the Rainbow sensors is specified in the following table.

Sensor	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
Rainbow DC-I	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3bpm
Rainbow DC-IP	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3bpm

SpCO accuracy is specified as ±3% under no motion conditions.

SpMet accuracy is specified as ±1% under no motion conditions.

Masimo SpO₂ Sensors

The Rad-57 may also use standard Masimo LNOP, LNOPv and LNCs SpO₂ sensors, when used with Red PC and LNCs Cabels respectively. This will allow the Rad-57 to work as a Masimo SET pulse oximeter without the carboxyhemoglobin or methemoglobin measurement.

Select the appropriate patient cable to attach the LNOP, LNOPv or LNCs sensor to the device.

RED REUSABLE SENSORS

(Red DCI sensors must be used in conjunction with Red PC cables)

SpO₂ and pulse rate accuracy for the Red sensors is specified in the following table.

Sensor	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
Red DC-I	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3bpm
Red DC-IP	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3bpm

LNOP® REUSABLE SENSORS

(LNOP sensors must be used in conjunction with Red PC cables)

Sensor	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNOP DC-I	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP DC-IP	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP YI	> 1 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	N/A	N/A
LNOP TC-I	> 30 kg	± 3.5%	N/A	± 3 bpm	N/A	± 3.5%	± 3 bpm
LNOP DC-195	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP DCSC	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP TF-I	> 30 kg	± 2%	N/A	± 3 bpm	N/A	± 2%	± 3 bpm

Note: The LNOP TF-I and TC-I sensors were not validated under motion conditions.

LNOP® ADHESIVE SENSORS

(LNOP sensors must be used in conjunction with Red PC cables)

Sensor	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNOP Adt	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP Adx	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP Pdt	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP Pdx	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP Neo	< 10 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNOP NeoPt	< 1 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNOP Neo-L	< 3 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNOP NeoPL	> 40 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNOP NeoPnL	< 1 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNOP Inf-L	3 - 20 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm

LNOP® SPECIALTY SENSORS

(LNOP sensors must be used in conjunction with Red PC cables)

Sensor	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNOP Blue	2.5 - 30 kg	60 - 80% ± 4%	N/A	± 3 bpm	N/A	± 3%	± 3 bpm
LNOP		70 - 100% ± 3.3%	N/A	± 3 bpm	N/A	± 3%	± 3 bpm
LNOP Newborn Neo	< 3 kg	80 - 100% ± 3%	N/A	± 3 bpm	N/A	± 3%	± 3 bpm
LNOP Newborn Neo	< 3 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNOP Newborn Inf	3 - 20 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm

LNCs™ REUSABLE SENSORS

(LNCs sensors must be used in conjunction with Red LNC cables)

SENSOR	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNCs DC-I	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs DC-IP	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs TC-I	> 30 kg	± 3.5%	N/A	± 3 bpm	N/A	± 3.5%	± 3 bpm
LNCs TF-I	> 30 kg	± 2%	N/A	± 3 bpm	N/A	± 2%	± 3 bpm

Note: The LNCs TF-I and TC-I sensors were not validated under motion conditions.

LNCs™ ADHESIVE SENSORS

(LNCs sensors must be used in conjunction with Red LNC cables)

SENSOR	Weight Range	Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
		No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNCs Adtx	> 30 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs Pdx	10 - 50 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs Inf-L	3-20 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs Neo-L	< 3 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm
LNCs NeoPnL	> 40 kg	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm
LNCs NeoPnL	< 1 kg	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm

LNOPV™ ADHESIVE SENSORS

(LNOPV sensors must be used in conjunction with Red PC cables)

SENSOR	Weight Range		Saturation Accuracy		Pulse Rate Accuracy		Low Perfusion Accuracy	
	< 3 kg	> 3 kg	No Motion	Motion	No Motion	Motion	Saturation	Pulse Rate
LNOPV Ne	± 3%	± 3%	± 3 bpm	± 5 bpm	± 3%	± 3 bpm	± 3%	± 3 bpm
LNOPV In	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm	± 2%	± 3 bpm
LNOPV Ad	± 2%	± 3%	± 3 bpm	± 5 bpm	± 2%	± 3 bpm	± 2%	± 3 bpm

SENSOR ACCURACY

Refer to Section 7, *Specifications for SpO₂, SpCO, SpMet and pulse rate accuracy*, unless otherwise specified in the tables above.

Accuracy specified when used with Masimo Rainbow SET technology pulse CO-Oximetry monitors or with licensed Masimo SET pulse oximetry modules during no motion. Numbers represent ± 1 standard deviation. Plus or minus one standard deviation represents 68% of the population. SpO₂ accuracy from 70% to 100%. Pulse Rate accuracy from 25 to 240bpm. Carboxyhemoglobin accuracy (SpCO) from 1 to 40%. Methemoglobin accuracy (SpMet) from 1 to 15%.

CLEANING AND REUSE OF MASIMO REUSABLE SENSORS AND CABLES

Reusable sensors and patient cables can be cleaned per the following procedure:

- Remove the sensor from the patient.
- Disconnect the sensor from the patient cable.
- Disconnect the patient cable from the monitor.
- Wipe the entire sensor and/or patient cable clean with a 70% isopropyl alcohol pad.
- Allow to air dry thoroughly before returning it to operation.

CAUTION: CAREFULLY ROUTE PATIENT CABLES TO REDUCE THE POSSIBILITY OF PATIENT ENTANGLEMENT OR STRANGULATION.

REATTACHMENT OF SINGLE USE ADHESIVE SENSORS

- Single use sensors may be reapplied to the same patient if the emitter and detector windows are clear and the adhesive still adheres to the skin.

NOTE: *If the sensor fails to track the pulse consistently, the sensors may be incorrectly positioned. Reposition the sensor or choose a different monitoring site.*

CAUTION: DO NOT REPROCESS ANY MASIMO SINGLE USE SENSORS.

Introduction

This chapter covers how to test the operation of the Rad-57, how to properly clean the Rad-57 Pulse CO-Oximeter, how to replace the batteries and how to obtain service.

Under normal operation, no internal adjustment or recalibration is required.

WARNING: BEFORE CLEANING THE OXIMETER, ALWAYS TURN IT OFF AND REMOVE THE BATTERIES.

Cleaning

To clean the display panel, use a cotton swab moistened with 70% isopropyl alcohol and gently wipe the panel.

To clean the outer surface of the oximeter, use a soft cloth dampened with a mild soap and water. Do not allow liquids to enter the interior of the instrument.

CAUTIONS:

- DO NOT AUTOCLAVE, PRESSURE STERILIZE, OR GAS STERILIZE THIS OXIMETER.
- DO NOT SOAK OR IMMERSE THE MONITOR IN ANY LIQUID.
- USE THE CLEANING SOLUTION SPARINGLY. EXCESSIVE SOLUTION CAN FLOW INTO THE MONITOR AND CAUSE DAMAGE TO INTERNAL COMPONENTS.
- DO NOT TOUCH, PRESS, OR RUB THE DISPLAY PANELS WITH ABRASIVE CLEANING COMPOUNDS, INSTRUMENTS, BRUSHES, ROUGH-SURFACE MATERIALS, OR BRING THEM INTO CONTACT WITH ANYTHING THAT COULD SCRATCH THE PANEL.
- DO NOT USE PETROLEUM-BASED OR ACETONE SOLUTIONS, OR OTHER HARSH SOLVENTS, TO CLEAN THE OXIMETER. THESE SUBSTANCES ATTACK THE DEVICE'S MATERIALS AND DEVICE FAILURE CAN RESULT.

Refer to Section 8, *Cleaning and Reuse of Masimo Reusable Sensors and Cables* for cleaning instructions of the sensor and patient cables.

BATTERY REPLACEMENT

The Rad-57 is powered by 4 "AA" alkaline batteries. Do not use any other type of batteries or power source to run the device. The battery compartment is accessed from the back of the device. To replace the batteries first remove the battery cover by depressing the small rectangular button at the bottom of the cover, and sliding the cover down off the bottom of the device. Remove the batteries and install new batteries in the directions indicated by the battery icons inside the battery compartment. Replace the battery cover by sliding it back up from the bottom of the device until the rectangular locking button snaps back into position.

WARNING: FAILURE TO REPLACE BATTERIES PROMPTLY AFTER A LOW BATTERY ALARM MAY RESULT IN THE CO-OXIMETER SHUTTING DOWN AND LEAVING THE PATIENT IN AN UNMONITORED CONDITION.

WARNING: USE ONLY ALKALINE BATTERIES. USE OF NON ALKALINE BATTERIES MAY AFFECT THE ACCURACY OF THE BATTERY LEVEL METER.

WARNING: USE OF BATTERIES WITH A CELL VOLTAGE OF MORE THAN 1.5V COULD CAUSE DAMAGE TO THE RAD-57.

WARNING: EFFECTIVE BATTERY LIFE WILL BE REDUCED WHEN OPERATING THE INSTRUMENT BELOW 5°F (-15° C) DUE TO ALKALINE BATTERY TECHNOLOGY.

Battery charge level is indicated by four LED indicators at the bottom of the front panel. All four indicators will be lit when the batteries are full, with fewer indicators being lit as the batteries lose their charge. When less than ten (10) percent battery life remains, the final battery indicator will begin to flash and an audible alarm will sound. Be sure to follow local regulations in regards to battery disposal.

Performance Verification

To test the performance of the Rad-57 Pulse CO-Oximeter following repairs or during routine maintenance, follow the procedure outlined in this section. If the Rad-57 fails any of the described tests, discontinue its use and correct the problem before returning the unit back to the user.

Before performing the following tests verify or install new batteries into the Rad-57 Handheld. Also disconnect any patient cables, pulse oximetry probes or serial cables from the instrument.

POWER-ON SELF-TEST

1. Turn the monitor on by depressing the Power Button. For about 5 seconds all available LEDs are illuminated and a brief beep tone sounds.
2. The oximeter begins normal operation.

KEY PRESS BUTTON TEST

1. With the exception of the Power Button, press each button and verify that the oximeter acknowledges each key-press with an audible beep tone or by indicating a change on the display.

ALARM LIMIT TEST

1. With the monitor turned on, select the Menu Access key and enter the Alarm menu. Change the High Saturation Alarm parameter to a value two points below the currently selected value, and accept the change.
 2. Verify that the newly set parameter is shown on the Saturation Alarm Limit Display, next to the SpO₂ or pulse rate measurement display.
 3. Return the High Saturation Alarm parameter to its original setting.
 4. Repeat steps 1 to 3 with the Low Saturation Alarm parameter.
 5. Repeat steps 1 to 3 with the High Pulse Rate Alarm parameter.
 6. Repeat steps 1 to 3 with the Low Pulse Rate Alarm parameter.
 7. Perform the above steps for the SpCO and SpMet alarm limits.
 8. Reset the alarm limits again to the original settings.

LED BRIGHTNESS

1. With the monitor turned on, select menu level 3 (see Section 4, *Setup Menu Level 3 - LED Brightness and Factory Defaults*) and use the Up and Down Arrow keys to cycle through all 4 brightness levels.
 2. Exit the Menu system by pressing the Mode/Enter key or waiting for the normal time-out.

Service and Repair

REPAIR POLICY

Masimo or an authorized Service Department must perform warranty repair and service. Do not use malfunctioning equipment. Have the unit repaired.

Please clean contaminated/dirty equipment before returning, following the cleaning procedure described in Section 9, *Cleaning*. Make sure it is fully dry before packing the equipment.

To return the Rad-57 unit for service, please follow the Return Procedure.

WARNING: DO NOT REMOVE THE COVER OF THE MONITOR EXCEPT FOR BATTERY REPLACEMENT. AN OPERATOR MAY ONLY PERFORM MAINTENANCE PROCEDURES SPECIFICALLY DESCRIBED IN THIS MANUAL. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL TRAINED IN THE REPAIR OF THIS EQUIPMENT.

RETURN PROCEDURE

Please clean contaminated/dirty equipment before returning and make sure it is fully dry before packing the equipment. Call Masimo at 800-326-4890 and ask for Technical Support. Ask for an RMA number. Package the equipment securely – in the original shipping container if possible – and enclose or include the following information and items:

- A letter describing in detail any difficulties experienced with the Pulse CO-Oximeter. Please include the RMA number in the letter.
- Warranty information – a copy of the invoice or other applicable documentation must be included.
- Purchase order number to cover repair if the oximeter is not under warranty, or for tracking purposes if it is.
- Ship-to and bill-to information.
- Person (name, telephone/telex/fax number, and country) to contact for any questions about the repairs.
- A certificate stating the oximeter has been decontaminated for bloodborne pathogens.

Return Rad-57 Pulse CO-Oximeter to the following shipping address:

For USA and Asia Pacific:	For Europe:
Masimo Corporation	Masimo Europe Limited
40 Parker	304 RN6, Le Bois des Cotes 2
Irvine, California 92618	69760 Limonest
Tel.: 949-297-7000	France
Fax.: 949-297-7001	

Warranty

Masimo warrants to the initial purchaser that each new Pulse CO-Oximeter will be free from defects in workmanship or materials for a period of one (1) year from the date of purchase. Masimo's sole obligation under this warranty is to repair or replace any product that Masimo deems to be covered under warranty with a repaired or a replacement Pulse CO-Oximeter.

Batteries are not warranted.

To request a replacement under warranty, contact Masimo for a returned goods authorization. If Masimo determines that a product must be replaced or repaired under warranty, it will be replaced or repaired and the cost of shipment covered. All other shipping costs shall be the responsibility of the purchaser.

Exclusions

This warranty does not extend to any product that has been subject to misuse, neglect or accident; that has been damaged by causes external to the Product; that has been used in violation of the operating instructions supplied with the product. The warranty does not extend to any product that has been connected to an unlicensed instrument system, modified accessories or any unit that has been disassembled or reassembled by anyone but an authorized Masimo agent.

THIS WARRANTY, TOGETHER WITH ANY OTHER EXPRESS WRITTEN WARRANTY THAT MAY BE ISSUED BY MASIMO IS THE SOLE AND EXCLUSIVE WARRANTY AS TO MASIMO'S PRODUCTS. THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY ORAL OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MASIMO SHALL NOT BE LIABLE FOR ANY INCIDENTAL, SPECIAL OR CONSEQUENTIAL LOSS, DAMAGE OR EXPENSE DIRECTLY OR INDIRECTLY ARISING FROM THE USE OR LOSS OF USE OF ANY PRODUCTS.

End-user license agreement

THIS DOCUMENT IS A LEGAL AGREEMENT BETWEEN YOU, THE "PURCHASER," AND MASIMO CORPORATION ("MASIMO"). IF YOU DO NOT AGREE TO THE TERMS OF THIS AGREEMENT, PROMPTLY RETURN THE ENTIRE PACKAGE, INCLUDING ALL ACCESSORIES, IN THEIR ORIGINAL PACKAGE, WITH YOUR SALES RECEIPT TO MASIMO FOR A FULL REFUND.

1. Grant of License: In consideration of payment of the license fee, which is part of the price paid for this product, MASIMO grants to Purchaser a nonexclusive, nontransferable license, without right to sublicense, to use the copy of the incorporated software/firmware and documentation in connection with Purchaser's use of the Masimo Products for their labeled purpose. MASIMO reserves all rights not expressly granted to Purchaser.

2. **Ownership of Software/Firmware:** Title to, ownership of, and all rights and interests in, any MASIMO software and/or firmware and the documentation, and all copies thereof, remain at all times vested in MASIMO Corporation, licensor to MASIMO, and they do not pass to Purchaser.
3. **Assignment:** Purchaser shall not assign or transfer this License, in whole or in part, by operation of law or otherwise, without MASIMO's prior written consent; any attempt without such consent, to assign any rights, duties or obligations arising hereunder shall be void.
4. **Copy Restrictions:** The software/firmware and the accompanying written materials are copyrighted. Unauthorized copying of the software, including software that has been modified, merged, or included with other software, or other written materials is expressly forbidden. You may be held legally responsible for any copyright infringement that is caused or incurred by your failure to abide by the terms of this license. Nothing in this license provides any rights beyond those provided by 17 U.S.C. §117.
5. **Use Restriction:** As the Purchaser, you may physically transfer the products from one location to another provided that the software/firmware is not copied. You may not electronically transfer the software/firmware from the products to any other device. You may not disclose, publish, translate, release or distribute copies of the software/firmware or accompanying written materials to others. You may not modify, adapt, translate, reverse engineer, decompile, disassemble, or create derivative works based on the software/firmware. You may not modify, adapt, translate, or create derivative works based on the written materials without the prior written consent of MASIMO.
6. **Transfer Restrictions:** The software/firmware is licensed to the Purchaser, and may not be transferred to anyone, except other end-users, without the prior written consent of MASIMO. In no event may you transfer, assign, rent, lease, sell, or otherwise dispose of the software/firmware or the products on a temporary basis.
7. **Beneficiary:** Masimo Corporation is a Beneficiary of this Agreement and has the right to enforce its provisions.
8. **U.S. Government Rights:** If you are acquiring software (including the related documentation) on behalf of any part of the United States Government, the following provisions apply: the software is deemed to be "commercial software" and "commercial computer software documentation," respectively pursuant to DFAR Section 227.7202 FAR 12.212, as applicable. Any use, modification, reproduction, release, performance, display or disclosure of the software (including the related documentation) by the U.S. Government or any of its agencies shall be governed solely by the terms of this Agreement and shall be prohibited except to the extent expressly permitted by the terms of this agreement.

Accessories

PART NUMBER	DESCRIPTION
1842	Rubber protective boot, Grey
1980	Rubber protective boot, Yellow
1981	Rubber protective boot, Red
1982	Rubber protective boot, Orange
2097	Rubber protective boot, Royal Blue
2098	Rubber protective boot, Light Blue
2099	Rubber protective boot, Pink
2208	Protective carrying case, black
2209	Protective carrying case, red
13158	Nylon protective carrying case
1908	CD, TrendCom Software
2063	PRONTO Cable, Rad-57, 20-pin



www.masimo.com

Instruments and sensors containing Masimo Rainbow SET technology are identified with the Masimo Rainbow SET logo.



© 2006 Masimo Corporation. Masimo, SET, LNCS, LNOP, Rainbow, SpCO, Signal IQ, Discrete Saturation Transform, DST and FastSat are registered trademarks of Masimo Corporation. Rad-57, Rad-57c, Rad-57m, Rad-57cm, SIQ, LNOPv, APOD, SpMet and Pulse CO-Oximeter are trademarks of Masimo Corporation.